
Now Introducing: All-inclusive Thermal Cycler Calibration Service

All-inclusive thermal cycler calibration service:

On-site PCR and qPCR thermal cycler calibration service


- No shipping of thermocycler
- Minimal down time

Comprehensive calibration certificate and report

- Overview charts and data summary
- Detailed charts and data per temperature step

Calibration data interpretation support

- Expert guidance how to interpret and incorporate calibration data



Accredited calibration service

CYCLERtest has been granted accreditation in 2009 by the Dutch Accreditation Council RvA. The accreditation is based on assessment against the requirements of the ISO/IEC 17025:2005. CYCLERtest has demonstrated to be technical competent, to be capable to generate technical valid results and work according to a management system. CYCLERtest is authorized to perform on-site PCR and qPCR thermal cycler calibrations under accreditation provided that the CYCLERtest continues to meet the requirements. The accreditation is registered under number K145 by the Dutch Accreditation Council

Accredited testing and medical laboratories

According to requirement 5.6.2.1 of ISO/IEC 17025:2005 and 5.6.3 of ISO/IEC 15189:2007 accredited testing or medical laboratories are required to calibrate their thermal cyclers traceable to the International System of Units (SI). Through accreditation CYCLERtest guarantees traceability to the SI.



Headquarters

CYCLERtest BV
Rotscherweg 61
6374 XW Landgraaf
The Netherlands

Phone: +31 45 533 8733
Fax: +31 45 533 8794
Website: www.cyclertest.com
Email: info@cyclertest.com